

trident

Phosphor Plates Advantages

- Flexible
- Cordless
- Safe
- Long Lasting
- Easy to positioning
- Reusable
- Thin



# READER

Phosphor Plate Scanner

# READER

Phosphor Plate Scanner

## PSP Scanner Trident Reader

### Simple

Connect Reader directly to your computer and start scanning; once you insert the plate into the slot, Reader automatically recognizes the plate's size and HD images will be in your computer in few seconds. Three LED lights will indicate the status of the process.

Reader has an inner power supply adaptor to work from 100 to 220 V.

### Sensor TDI CCD

Reader does not need big spaces for installation. The whole system is wonderfully compact, measuring 17 x 19 x 30 cm only.


### Fast

Reader delivers highly precise digital images to any PC in less than 20 seconds: Only ten seconds to scan and five seconds to clean the plates for a new image.

### Affordable

Reader is a cost-effective option for any practice. It perfectly works with your traditional AC or DC X-Ray generator, providing high quality images for accurate diagnosis.

Reader enables getting clear HD images repeatedly using the same plate, with no film processing extra cost: No chemicals, film processor or dark room needed.



Sizes 0, 1, 2 and 3



## xrddental

Dental Equipment - [www.xrddental.com](http://www.xrddental.com)

XRD Dental Equipment - Orlando FL 32836  
+ 1 (407) 319-1348 - [info@xrddental.com](mailto:info@xrddental.com)